

# Reliable Stainless Steel Products

- √ Cabinets
- √ Tables / Carts
- √ Bassinets
- √ Mayo Instrument Stands
- √ I.V. Stands
- √ Receptacles
- √ Stools



*UMF Quality, Design and Innovation*

**EXPECT MORE**

# STAINLESS STEEL STORAGE and SUPPLY CABINETS

---



**SS7816**  
Storage and Supply  
Cabinet

## SS7816 • Storage and Supply Cabinet

Five adjustable stainless steel shelves, two sliding glass-paneled doors, leveling glides

Option: SS7816-503 Lock (Factory installed only)

Overall Dimensions: 47" W x 84" H x 16" D

## SS7840 • Storage and Supply Cabinet

Five adjustable stainless steel shelves, two hinged glass-paneled doors, T-handle lock, leveling glides

Overall Dimensions: 36" W x 77" H x 16" D

## SS7834 • Instrument Cabinet

Two hinged glass-paneled doors, T-handle lock, five adjustable glass shelves, leveling glides

Option: SS7834-380 Stainless steel shelves

Overall Dimensions: 36" W x 78" H x 16" D

*All stainless steel storage and supply cabinets are available with solid doors. Please call UMF Customer Service for details.*



**SS7834**  
Instrument Cabinet



**SS7840**  
Storage and  
Supply Cabinet



*UMF Quality, Design and Innovation*

**EXPECT MORE**

[www.umf-exam.com](http://www.umf-exam.com)

# NARCOTICS STORAGE CABINETS

---

## SS7780 • Double Door Narcotics Storage Cabinet

Outer door is 20 gauge, inner door is 14 gauge, separate locks, different keys, three adjustable shelves

Option: SS7780-530 Locker warning light (factory installed only)

Overall Dimensions: 18" W x 24" H x 12" D



**SS7780**

Double Door  
Narcotics Storage Cabinet

## SS7781 • Double Door Narcotics Storage Cabinet

Outer door is 20 gauge, inner door is 14 gauge, separate locks, different keys, two adjustable shelves

Option: SS7781-530 Locker warning light (factory installed only)

Overall Dimensions: 12" W x 15" H x 8" D



**SS7781**

Double Door  
Narcotics Storage Cabinet

## SS7782 • Single Door / Single Lock Narcotics Storage Cabinet

Single door, can be wall mounted or installed into right-hand corner of wall cabinet, extra wide end permits door to swing open fully, one adjustable shelf

Option: SS7782-530 Locker warning light (factory installed only)

Overall Dimensions: 14" W x 12" H x 8" D

## SS7784 • Single Door / Single Lock Narcotics Storage Cabinet

Single door, single lock, no shelf

Option: SS7784-530 Locker warning light (factory installed only)

Overall Dimensions: 12" W x 7.75" H x 7" D



**SS7782**

Single Door / Single Lock  
Narcotics Storage Cabinet



**SS7784**

Single Door / Single Lock  
Narcotics Storage Cabinet

# DRESSING CART

---

## SS8214 • Dressing Cart and Utility Table (Knock Down)

Knock down construction, guard rails are stainless steel tubing .5" diameter. Basin and pail ring supports are 1" x .25" stainless steel with rubber bumpers to prevent rattling. 5" rubber tire casters, two swivel and two rigid, with bumper guards. Standard equipment includes 7.5 quart basin and 13 quart pail, all stainless steel.

Overall Dimensions: 22.25" W x 52.125" L x 39.75" H

Shelf Dimensions: 21.75" W x 37.5" L



**SS8214**  
Dressing Cart and  
Utility Table

## UTILITY TABLES (KNOCK-DOWN)

---

Table supported by four uprights and four-sided guard rails. Table top and shelf edges are bent and hemmed. 3" ball-bearing swivel casters.

### SS8090 • Utility Table - Knock Down

Overall Dimensions: 33" W x 34" H x 18" D

### SS8096 • Utility Table - Knock Down

Overall Dimensions: 20" W x 34" H x 16" D



**SS8090**  
Utility Table - Knock Down



**SS8096**  
Utility Table - Knock Down

# UTILITY TABLES

All stainless steel construction. Three-sided guard rails, 2" ball-bearing swivel casters, all edges bent and hemmed.

## SS8120 • Utility Table (Two drawers, one shelf)

Overall Dimensions: 32" W x 32" H x 16" D



SS8120  
Utility Table

## SS8150 • Utility Table (Four drawers)

Overall Dimensions: 20" W x 35.25" H x 16" D



SS8150  
Utility Table -  
4 Drawers

## SS8151 • Utility Table (Three drawers, one shelf)

Overall Dimensions: 20" W x 35.25" H x 16" D

## SS8152 • Utility Table (Two drawers, one shelf)

Overall Dimensions: 20" W x 35.25" H x 16" D

## SS8153 • Utility Table (One drawer, one shelf)

Overall Dimensions: 20" W x 35.25" H x 16" D

### Available Options:

- 505 Lock for single drawer (Factory installed only)
- 480 rubber crutch tips (Height will be 2 inches lower)



SS8151  
Utility Table - 3 Drawers



SS8152  
Utility Table - 2 Drawers



SS8153  
Utility Table -1 Drawer

# INSTRUMENT / BACK TABLES with SHELF OR H-BRACE

Made with 20 gauge stainless steel. Tops and shelves are hemmed and formed for extra perimeter strength. Stiffeners are included on larger size tables. All tables are 34" in height with 2" casters. 19" space between shelves.

with shelves	Legs = 16 Gauge SS Tubing		Top w/ Stiffener	Shelf w/ Stiffener
SS8000	24 x 72	1.25" O.D.	x	x
SS8004	24 x 60	1.25" O.D.	x	x
SS8008	24 x 48	1" O.D.	x	x
SS8010	20 x 36	1" O.D.	x	x
SS8012	18 x 33	1" O.D.	x	x
SS8014	20 x 48	1" O.D.	x	x
SS8016	16 x 30	1" O.D.	x	x
SS8018	16 x 20	1" O.D.	x	x

with h-brace	Legs = 16 Gauge SS Tubing		Top w/ Stiffener	Shelf w/ Stiffener
SS8030	24 x 72	1.25" O.D.	x	
SS8038	24 x 48	1.25" O.D.	x	
SS8042	20 x 36	1" O.D.	x	

Options: -480 rubber tips (32" H); .375 diameter guard rail (contact UMF for pricing)



**SS8010**  
Instrument / Back Table  
20x36



**SS8012**  
Instrument / Back Table  
18x33



**SS8014**  
Instrument / Back Table  
20x48



**SS8018**  
Instrument / Back Table  
16x20

# BASSINETS

Basket supports with insert bends allow basket to be moved to horizontal or elevated positions. Mounted on 3" casters, tubular framework is 1" O.D. x 16 gauge stainless steel, braced with .5" round rods. Cabinets are 20 gauge stainless steel. Guard rails are .375". Baskets and mattresses are optional on all models.

## SS8525 • Bassinet

One shelf, one drawer  
Overall Dimensions: 32" W x 36.25" H x 17" D

## SS8528 • Bassinet

Two doors  
Overall Dimensions: 32" W x 36.25" H x 17" D

## SS8538 • Bassinet

One shelf  
Overall Dimensions: 32" W x 36.25" H x 17" D

## 8545 • Basket

Plastic basket, one piece construction, crystal clear for maximum visibility, name card holder  
Overall Dimensions: 16.25" W x 9.75" H x 30" L

## 8546 • Mattress

Foam, covered with vinyl sheeting  
Overall Dimensions: 13" W x 26" L x 1.5" thick



\*Shown with 8545 Basket and 8546 Mattress

## MAYO INSTRUMENT STANDS



Non-tip weighted flat base slides easily under O.R. table. Removable stainless steel tray is standard equipment. Mounted on 2" ball-bearing swivel casters.

Model	Height	Tray	Adjustment
SS8310	39.5 - 58.5"	13 x 19	Foot pedal
SS8311	39.5 - 58.5"	16 x 21	Foot pedal
SS8318	33 - 58"	13 x 19	Hand wheel
SS8319	33 - 58"	16 x 21	Hand wheel



## I.V. STANDS

24" diameter base, sturdy, functional design mounted on five legged, non-tip aluminum base for stability. Fitted with 2" swivel casters. Telescopic, stainless steel rod adjusts by positive locking hand wheel. 1" diameter tubing.



Model	Hook	Extended Height	Adjustment
SS8342	4	92.5"	Hand wheel
SS8343	2	92.5"	Hand wheel

# BASIN STANDS

---

Tubular framework fitted with basin holder, 13" stainless steel basin (7 quart). Rubber bumper inside basin holder eliminates rattling and scarring of basin. Mounted on 2" ball-bearing swivel casters. Stands 32" high. H-Brace or shelf for rigid construction.

SS8360 • Double Basin Stand with Shelf

SS8361 • Double Basin Stand with H-Brace

SS8365 • Single Basin Stand with Shelf

SS8366 • Single Basin Stand with H-Brace

SS8366  
Single Basin Stand  
with H-Brace

SS8360  
Double Basin Stand  
with Shelf



# KICK BUCKET, SPONGE RECEPTACLE, WASTE RECEPTACLE

---

SS8350  
Kick Bucket



All-welded construction, non-tip, solid balance, heavy duty rubber bumper, rubber cushions inside frame prevents pail from rattling.

SS8350 • Kick Bucket

13 quart stainless steel pail mounted on stainless steel frame with four 2" ball-bearing swivel casters

SS8355 • Sponge Receptacle

10 quart basin mounted on stainless steel frame with four 2" ball-bearing swivel casters



SS8355  
Sponge Receptacle

SS6849 • Waste Receptacle

Foot-operated lid, 20 quart capacity, stainless steel inner pail with handle

Overall Dimensions: 13" W x 18" H x 12" D

SS6849  
Waste Receptacle



# STEP STOOLS

---

All welded construction. Non-tip offset legs, 16 gauge x 1" O.D. SS tubing, corrugated rubber matting cemented to platform. Platform has continuous 1.25" x .125" stainless steel band. Rubber floor tips.

## SS8370 • Foot Stool ~ 2 steps

Overall Dimensions: 18" W x 15.625" H x 25" D



SS8370  
Foot Stool ~ 2 steps

## SS8372 • Foot Stool ~ 1 step

Overall Dimensions: 30" W x 7.75" H x 12" D



SS8372  
Foot Stool ~ 1 step

## SS8374 • Foot Stool ~ 1 step

Overall Dimensions: 18" W x 7.75" H x 12" D

## SS8376 • Foot Stool ~ 1 step

Overall Dimensions: 14" W x 7.75" H x 10" D

## SS8378 • Foot Stool w/ Handrail ~ 1 step (not shown)

Overall Dimensions: 18" W x 7.75" H x 12" D Handrail 35.75"

## SS8380 • Stackable Foot Stool

Overall Dimensions: 18" W x 6" H x 12" D

## SS8381 • Transport Cart

Overall Dimensions: 20.5" W x 8.5" H x 15" D Handrail 36"



SS8380  
Stackable Foot Stool

SS8381  
Transport Cart

SS8374  
Foot Stool ~ 1 step

# REVOLVING and FIXED SEAT STOOLS

---

Non-slip material mounted into recess of stool seat, flush finish. Seat diameter is 14", seat adjusts with fast-action 1" O.D. seat adjustment screw. Welded footrest or cross bar with rubber tips.

## SS6700 • Revolving Seat Stool

with backrest, adjusts from 25" to 31", ring foot rest is 16.25" square

Option: SS6700-112 • 2" ball-bearing swivel caster (factory-installed only)

## SS6701 • Revolving Seat Stool (not shown)

without backrest, adjusts from 25" to 31", ring foot rest is 16.25" square

Option: SS6701-112 • 2" ball-bearing swivel caster (factory-installed only), -uph • upholstered seat cushion

## SS6704 • Revolving Seat Stool

adjusts from 19" to 25", ring foot rest

Option: -uph • upholstered seat cushion

-112 casters

## SS6725 • Fixed Seat Stool

cross brace, rubber tips only, height is 18", seat diameter is 14"

Option: SS6725-112 • 2" ball-bearing swivel caster (factory-installed only)



**SS6700**

Revolving Seat Stool



**SS6704**

Revolving Seat Stool  
Shown with optional  
-112 casters



**SS6725**

Fixed Seat Stool



UNITED METAL FABRICATORS

**United Metal Fabricators, Inc.**  
1316 Eisenhower Blvd.  
Johnstown, PA 15904  
800-638-5322 • 814-266-1870 fax  
[www.umf-exam.com](http://www.umf-exam.com)