

Unsurpassed quality and functionality

- UMF Ultra Comfort Seamless Top ■ Easy Patient Positioning
- Large Storage Drawer ■ 600 Pound Capacity
- 19" Low Access ■ 3 Year Limited Warranty



4040



4070

The 4040 and 4070 lower to a height of 19 inches from the floor and can support a patient weighing up to 600 pounds. The Ultra Comfort Seamless Top provides security and comfort for patients and also allows for easy cleaning and helps prevent cross-contamination.

Both models have a power hi-low feature which positions patients at the press of a button, and the 4040 also includes a fully articulating, power assisted back section allowing for simple and easy patient position changes.

An elevated leg rest pad is standard on both models, providing additional patient comfort.

Operating on heavy-duty roller bearing drawer glides, our drawers provide dependable, convenient operation for as long as you own your table.

Standard Features:

- UMF Ultra Comfort Seamless Top
- 13 Upholstery color options (see page 12)
- Removable stainless steel drain pan (3.5 qt. / 3.3 L)
- Stainless steel leg pullout
- Removable upholstered leg cushion
- Leg section elevates 10 degrees for patient comfort
- One touch, easy-glide stirrups
- Electro-mechanical actuator system
- Concealed paper roll holder (18"-21") (45.7cm-53.3cm)
- Front storage drawer
- Two function foot control (4040)
- One function foot control (4070)
- Easy release manual back section (4070)
- 600 lb. (272.2 kg.) patient weight capacity
- Pelvic tilt

Table Specifications:

- Height: min. 19" (48.3cm) at seat max. 31.5" (80cm) at seat
- Width: 27" (68.6cm)
- Length: 58.5" (148.6cm) 70" (177.8cm) with foot section extended
- Leg Rest: 20" (50.8cm) W x 12" (30.5cm) H
- Drawer: 17" (43.2cm) W x 13" (33cm) L x 4 3/4" (12.1cm) H
- Electrical: 115 VAC, 60 Hz
- Stirrups: 3 lateral positions

Options:

(see page 11 for popular option details)

- 011 Leg section booster pad (4040/4070)
- 035 Armboard assembly (contains one each 036, 038, 039)
- 036 Sliding bracket and lug (reversible) (required for 038)
- 038 Armboard (requires 039 accessory rails)
- 039 Accessory rails (factory installed)
- 063 Two function hand control (4040)
- 065 Vision screen
- 066 One function hand control (4070)
- 143 Easy-Glide locking casters
- 251 Bierhoff knee crutches
- 375 Safety rails (pair) (requires 2 ea. 039 accessory rails)
- 810 Restraint straps (each strap requires 2 ea. 039 accessory rails)
- 822 Stainless steel drain pan with spout
- 850 Drawer warmer 4040/4070 only (factory installed)
- 852 Hospital grade receptacle
- 9999 T-cable for simultaneous use of both hand and foot controls
- Fire retardant foam
- 220 VAC, 50/60 HZ