


Providing the capacity to lift 800 pounds

- UMF Ultra Comfort Seamless Top
- 19" Access
- 800 Pound Capacity
- Ultimate Stability
- 3 Year Limited Warranty



5005

- Power Hi-Low
- Power Back

The 5005 Bariatric Treatment Table provides the features you need, along with the comfort and stability your patients expect. The 5005 is the only bariatric power treatment table to feature a 30-inch wide and 70-inch long vacuum-formed UMF Ultra Comfort Seamless Top extendable to 80" with the articulating leg extension. The design allows for easy cleaning and helps prevent cross-contamination; while the thick, high density pad with built-in pillow provides the ultimate in patient comfort.

The 5005 can easily and safely lift up to 800 pounds from its 19-inch low patient access to its highest position at 35 inches. The low access makes this table easily accessible for wheelchair transfers, elderly patients and patients with special needs. Assisting patients onto this table is less stressful on medical staff (helps prevent back injuries), and the range from 19 to 34 inches allows the table to be adjusted to the doctor's ideal working height.

The 5005 has a fully articulating, power assisted back section as well as an articulating leg section that provides added comfort and support for taller patients.

Standard Features:

- UMF Ultra Comfort Seamless Top
- 13 Upholstery color options (see page 12)
- Built-in pillow
- 2 Function hand control
- Electro-mechanical actuator system
- Concealed paper roll holder (18" -21") (45.7cm-53.3cm)
- Large storage drawer (easily changed from patient left to patient right)
- Fully articulating 10" (25.4cm) leg section
- Accessory rails (at chest)
- 800 lb. (362.9 kg.) patient weight capacity

Table Specifications:

- Height: min. 19" (48.3cm) max. 35" (88.9cm)
- Width: 30" (76.2cm)
- Length: 70" to 80" (177.8cm to 203.2cm)
- Storage Drawer: 16" D x 19" W x 4.5" H (40.6cm x 48.3cm x 11.4cm)
- Electrical: 115 VAC, 60 HZ

Options:

- (see page 11 for popular option details)
- 036 Sliding bracket and lug (reversible) (required for 038)
 - 038 Armboard only (used with 036)
 - 039 Accessory rails (factory installed) (pair for seat only)
 - 062 Two function foot control
 - 375 Safety rails (pair) (requires 2 ea. 039 accessory rails)
 - 800 Paper cutter band
 - 810 Restraint straps
 - 9999 T-cable for simultaneous use of both hand and foot controls
 - Fire retardant foam
 - 220 VAC, 50/60 HZ