



5020

Our most economical power procedure table

- UMF Ultra Comfort Seamless Top
- 24" Access
- 400 Pound Capacity
- 3 Year Limited Warranty

The 5020 lowers to a height of 24 inches in a reassuring upright chair position.

The one-touch, four function hand control allows for precise, smooth operation. A four function programmable foot control is available that makes positioning the table in your favorite position quick and effortless.

The UMF Ultra Comfort Seamless Top provides security and comfort for patients and also allows for easy cleaning and helps prevent cross-contamination.

The 5020 has a capacity of 400 pounds and provides the versatility needed for OB-GYN, dermatology, proctology, minor surgery and urology procedures and exams.

Standard Features:

- UMF Ultra Comfort Seamless Top
- 13 Upholstery color options (see page 12)
- Articulating headrest
- Stainless steel leg pullout
- Removable stainless steel drain pan (1.9 qt. / 1.8 L)
- Two position upholstered kneeler/footrest
- One touch, easy-glide stirrups
- Electro-mechanical actuator system
- Concealed paper roll holder (18" -21") (45.7cm-53.3cm)
- Trendelenburg position
- Four function hand control (hi-low, back, leg rest, seat tilt)
- 400 lb. (181.4 kg.) patient weight capacity

Table Specifications:

- Height: min. 24" (61cm) at seat max. 36" (91.4cm) at seat
- Width: 26.5" (67.3cm)
- Length: 71.5" (181.6cm) with foot section extended 85" (215.9cm) with foot section and head rest extended
- Leg Rest: 16.25" (41.3cm) W x 9" (22.9cm) H
- Head Rest: 15" (38.1cm) W x 8.5" (21.6cm) H
- Electrical: 115 VAC, 60 Hz
- Stirrups: 3 lateral positions

Options:

(see page 11 for popular option details)

- 035 Armboard assembly (contains one each 036, 038, 039)
- 036 Sliding bracket and lug (reversible) (required for 038)
- 038 Armboard (requires 039 accessory rails)
- 039 Accessory rails (factory installed)
- 060 Four function programmable foot control
- 065 Vision screen
- 143 Easy-Glide locking casters
- 251 Bierhoff knee crutches
- 375 Safety rails (pair) (requires 2 ea. 039 accessory rails)
- 810 Restraint straps (each strap requires 2 ea. 039 accessory rails)
- 822 Stainless steel drain pan with spout
- 9999 T-cable for simultaneous use of both hand and foot controls
- Fire retardant foam
- 220 VAC, 50/60 HZ