

Our most economical power exam table

- UMF Ultra Comfort Seamless Top
- 400 Pound Capacity
- 3 Year Limited Warranty

The 5060 and 5080 both offer multiple positions for various exams and procedures. The 5080 is easily raised and lowered with its foot pedal control, while the 5060 utilizes a two function hand control for both hi-low and power back functions. Both models can also be used as a chair in their full upright positions.

The 5060 and 5080 feature a stainless steel leg section with a removable upholstered leg cushion. The leg section elevates 10 degrees for added patient comfort.

As with all of our power tables and chairs, the 5060 and 5080 utilize the vacuum-formed UMF Ultra Comfort Seamless Top which provides security and comfort for patients, allows for easy cleaning and helps prevent cross-contamination.



5060

- Power Hi-Low and Power Back
- Hand Held Control



5080

- Power Hi-Low
- Foot Operated Control

Standard Features:

- UMF Ultra Comfort Seamless Top
- 13 Upholstery color options (see page 12)
- Removable stainless steel drain pan (3.5 qt. / 3.3 L)
- Stainless steel leg pullout
- Removable upholstered leg cushion
- Leg section elevates 10 degrees for patient comfort
- One touch easy-glide stirrups
- Electro-mechanical actuator system
- Concealed paper roll holder (18" -21") (45.7cm-53.3cm)
- 2 Function hand control (5060)
- 1 Function foot control (5080)
- Easy release manual back section (5080)
- 400 lb. (181.4 kg.) patient weight capacity
- Pelvic tilt

Table Specifications:

- Height: min. 24" (61cm) at seat max. 36" (91.4cm) at seat
- Width: 26.5" (67.3cm)
- Length: 54.5" to 74" (138.4cm to 188cm)
- Electrical: 115 VAC, 60 HZ
- Stirrups: 3 lateral positions

Options:

(see page 11 for popular option details)

- 011 Leg section booster pad
- 035 Armboard assembly (contains one each 036, 038, 039)
- 036 Sliding bracket and lug (reversible) (required for 038)
- 038 Armboard (requires 039 accessory rails)
- 039 Accessory rails (factory installed)
- 062 2 function foot control (5060)
- 065 Vision screen
- 066 1 function hand control (5080)
- 143 Easy-Glide locking casters
- 251 Bierhoff knee crutches
- 375 Safety rails (pair) (requires 2 ea. 039 accessory rails)
- 810 Restraint straps (each strap requires 2 ea. 039 accessory rails)
- 822 Stainless steel drain pan with spout
- 845 Accessory drawer
- 9999 T-cable for simultaneous use of both hand and foot controls
- Fire retardant foam
- 220 VAC, 50/60 HZ