


neonatal scale



model 610-10-2



RICE LAKE[®]
HEALTH SOLUTIONS

To be the best by every measure[®]

www.ricelake.com/health

neonatal scale

THE 610-10 NEONATAL SCALE provides the unsurpassed ability to weigh an infant down to the gram. Combine this high-resolution accuracy with movement compensation and reweigh function and you have a winning solution.

The movement compensation feature of the Rice Lake Neonatal Scale improves accuracy and speed when weighing delicate newborns. Advanced speed means fast weighing, allowing the baby's quick return to the incubator.

Durable, easy to use and feature-rich, the Rice Lake Neonatal Scale gives healthcare professionals all the tools they need to provide quality infant care and precision measurement.

Capacity	0-5 lb (0-2 kg)	5-12 lb (2-6 kg)	12-30 lb (6-15 kg)
Graduations	0.05 oz (1 g)	0.1 oz (2 g)	0.2 oz (5 g)



Optional stainless steel cart.



Pictured with optional four-sided tray.

610-10-2 Standard Features



- Tray dimensions: 12 in W x 22 in L x 4¼ in H
- Removeable tray with built-in measuring tape (up to 21¼ in x 1/8 in or 55 cm x ¼ cm)
- 1 inch LCD display
- Buttons: reweigh, lb/kg, zero, on/off
- Reweigh function
- Low-battery indication
- Motiontrap (movement compensation technology)

Options/Accessories

- Four-sided clear tray with built-in measuring tape
- Stainless steel cart

Warranty

Two-year limited warranty

Your Rice Lake Weighing Systems distributor is:

RICE LAKE
HEALTH SOLUTIONS

To be the best by every measure®

230 W. Coleman St. • Rice Lake, WI 54868 • USA
TEL: 715-234-9171 • FAX: 715-234-6967 • www.ricelake.com

An ISO 9001 Registered Company © 2010 Rice Lake Weighing Systems PN 106420 4/10
Specifications subject to change without notice