



NOVATECH
MEDICAL

Procedure & Exam Lighting

LED NOVA 4

Precision Lighting. Exceptional Service.

NOVATECH
MEDICAL

Novatech Medical Inc.
H-3955 Isabelle Street,
Brossard, QC J4Y 2R2

Call : (450) 449-3608

Visit : www.novatechmedical.com

Email: info@novatechmedical.com



Focused illumination for precision and comfort

The **LED NOVA 4** is designed to deliver clear, consistent lighting for examinations and minor procedures. With high colour rendering ($\text{CRI} \geq 95$) and uniform light distribution, it ensures accurate tissue assessment and improved patient outcomes in any clinical setting.

Its compact, lightweight form and flexible arm system allow for easy positioning and reduced shadows—whether wall-mounted, ceiling-suspended, or mobile. The **LED NOVA 4** adapts seamlessly to your workflow, enhancing comfort for both patients and practitioners.

Smart design. Reliable performance.

Built with energy-efficient, long-lasting LED technology, the **LED NOVA 4** offers cool-touch operation and virtually maintenance-free performance. Backed by Novatech Medical's commitment to practical, ergonomic innovation, it's a dependable choice for any modern healthcare environment.

Better visibility. Better care. Every time.



TECHNICAL SPECIFICATIONS	LED NOVA 4
Light head diameter	60 cm
Number of LED's	39 LEDs
Light intensity at 1m	130 000 lux
Light field diameter	10-25 cm
Colour temperature	3500 ; 4000 ; 4500 ; 5000 ; 5500 K
Colour rendering index (Ra)	95
Colour rendering index (R9)	93
Lifespan of LED's	50 000 hours
Control	Light head control ; Sterile handle ; Optional wall control



Easily accessible ENDO & R9 (red light) mode



Intuitive COLOUR TEMPERATURE & LIGHT INTENSITY controls located directly on the light head



FOCUSABLE light field control using central handle & removable STERILIZABLE handle



Mounting options include CEILING, WALL, and ROLL-CART



**NOVATECH
MEDICAL**

Novatech Medical LED Nova 4 Procedure & Examination Light



NOVATECH MEDICAL

Learn more about Novatech Medical products at
www.novatechmedical.com or call us at **(450) 449-3608**