

SS6704

SS6700 Series Stools



Stainless Steel Revolving Stool with Foot Ring



*Stool shown with option 112 – Swivel Casters

*Mild corrosion on the stem of the stool can take place if subject to moisture (SS6704 only)

Standard Features

- High quality, welded stainless steel construction
- Welded foot rest
- Non-slip material mounted into recess of stool seat (flush finish)
- Non-slip rubber feet on stool legs (device shown with optional casters)
- Heavy duty threaded stem for seat height adjustment
- 3 year limited warranty

Optional Accessories

- **112** – Swivel Casters 2" (set of 4)
- **117** – Swivel Casters Locking 2" (set of 4)
- **SS6704NC** – Non-corrosive steel

Specifications

- **Seat diameter:** 14" (35.6cm)
- **Foot rest dimensions:** 16.25"x16.25" (41.3cm x 41.3cm)
- **Min. seat height:** 19" (48.3cm)
- **Max. seat height:** 25" (63.5cm)

