

SS8372

SS8300 Series Foot Stools



Stainless Steel One Step Foot Stool



Standard Features

- High quality, welded stainless steel construction
- 16 gauge, 1" (2.54cm) outside diameter stainless steel tubing
- 1.25" x .125" (3.2cm x .3cm) stainless steel band below step platform for rigidity
- Non-slip rubber feet on stool legs
- Non-slip rubber step treads
- 3 year limited warranty

Specifications

- **Width:** 30" (76.2cm)
- **Height:** 7.75" (19.7cm)
- **Depth:** 12" (30.5cm)

