

SS8525

SS8500 Series Bassinets



Bassinet, One Shelf, One Drawer with Basket and Mattress



Standard Features

- Heavy duty tubular stainless steel framework
- Full width stainless steel drawer
- Convenient storage shelf
- 3" (7.6cm) ball bearing swivel casters
- 3 year limited warranty
- Includes basket and mattress pad
- Name card holder

Optional Accessories

- 132 – Swivel casters locking 3" (set of 4)

Specifications

- **Height:** 35.375" (89.9cm)
- **Width:** 31.25" (79.4cm)
- **Depth:** 17" (43.2cm)

