

4040-650-200

Power Hi-Lo & Power Back Exam Tables



Standard Features

- Power Hi-Lo and power backrest
- 650lb (294.84kg) patient weight capacity
- Removable upholstered top with seamless antimicrobial guard
- 3 year limited warranty
- Low access ADA compliant 18.5" (±0.5") ((46.9cm)(±1.27cm))
- Two function foot control with Hi-Lo and power back
- Large front storage drawer
- Side storage drawer - reversible
- Dual electro-mechanical actuator system for efficient speed
- 12 gauge steel extension leg pull-out to accommodate all patient heights
- Seamless upholstered leg cushion elevates 10° for patient comfort
- Paper roll holder (18" – 21") (45.7cm – 53.3cm)
- One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment
- Pelvic tilt (7°)
- Magnetic hand control with OneTouch Patient Positioning™ — two programmable table positioning buttons, plus return-to-home
- Hospital grade receptacle 110V with child safety cover (patient left)
- Stainless steel drain pan (1.5" deep) (3.81cm)

Optional Accessories

- **014** – Leg extension booster pad to elevate leg section (*field installable*)
- **146** – OneTouch Wheelbase™ for infection control, ADA compliant
- **251** – Bierhoff knee crutches (*pair/field installable*)
- **469** – Welch Allyn Light bracket capabilities (NOT actual bracket) (*specify patient right/left*)
- **822** – Stainless steel drain pan with spout (*factory installed only*)
- **870** – Drawer warmer 110V (*factory installed only*)
- **4040-650-200F** – Fire retardant foam (CAL133)

Table Specifications

- **Height:**
Min. 18.5" (± .5") (46.9cm) (±1.27cm) at seat
Max. 37" at seat (93.9cm)
- **Width:** 27" (68.6cm)
- **Length:** 57" – 76" with foot section extended (144.7cm - 193cm)
- **Leg rest:** 21" W x 13" D (53.3cm W x 33cm D)
- **Front Storage Drawer:** 19" W x 12" D (48.26cm W x 30.48cm D)
- **Rear Storage Drawer:** 12.5" W x 24.5" D (31.75cm W x 62.23cm D)
- **Electrical:** 120 VAC, 60HZ
- **Stirrups:** Four lateral positions, Infinitely adjustable

Colors

