

4010

4000 Series Power Procedure / Exam Tables



Four Function Power Procedure Table



Standard Features

- Ultra comfort seamless top with antimicrobial guard
- ADA compliant
- 19" (48.3cm) low access seat height
- Power Hi-Lo
- Power back
- Power tilt
- Power leg section
- Four function hand control (hi-lo, back, leg rest, and seat tilt)
- Articulating headrest
- Stainless steel leg pullout
- Removable stainless steel drain pan (3.5qt./3.3L)
- Two position upholstered kneeler/footrest
- One touch, easy-glide stirrups with 3 lateral positions and friction lock length adjustment
- Concealed paper roll holder (18"-21") (45.7cm-53.3cm)
- Trendelenburg position
- Electro-mechanical actuator system
- 600lb (272.2kg) patient weight capacity
- 3 year limited warranty

Optional Accessories

- 035 – Arm board assembly (*specify patient right or left/factory installed*)
- 041 – Patient assist armrests (*pair/factory installed*)
- 045 – Accessory rail (*pair/specify head or seat section/factory installed*)
- 065 – Vision screen
- 143 – Easy-Glide locking casters (*set of four/ factory installed*)
- 251 – Bierhoff knee crutches (*pair*)
- 375 – Safety rails full assembly (*pair/factory installed*)
- 411 – Four function programmable foot control (*simultaneous use of foot and hand controls*)
- 810 – Restraint straps assembly (*pair/factory installed*)
- 822 – Stainless steel drain pan with spout (*factory installed*)
- 852 – Hospital grade receptacle
- 4010F – Fire retardant foam

Table Specifications

- **Height:** Min. 19" (48.3cm) at seat
Max.35" (88.9cm) at seat
- **Width:** 27" (68.6cm)
- **Length:** 71.5" (181.6cm) with foot section extended; 85" (215.9cm) with foot section and headrest extended
- **Leg rest:** 16.25" (41.3cm) W x 9" (22.9cm) L
- **Head rest:** 15" (38.1cm) W x 8.5" (21.6cm) H
- **Electrical:** 120 VAC, 60 HZ or 220/240VAC, 60HZ
- **Stirrups:** Infinitely adjustable

Colors

