

5005

5000 Series Bariatric Treatment Tables



Two Function Bariatric Power Table



Standard Features

- Ultra comfort seamless top with antimicrobial guard and built in pillow
- ADA Compliant with 19" (48.3cm) low access
- All-steel heavy duty construction
- Power Hi-Lo
- Two function hand control
- Concealed paper roll holder (18" – 21") (45.7cm – 53.3cm)
- Large, full extension storage drawer on roller bearing drawer guides
- 10" (25.4cm) leg section
- Electro-mechanical actuator system
- 800lb (363kg) weight capacity
- 3 year limited warranty

Optional Accessories

- **035** – Armboard assembly (*specify patient right or left/factory installed*)
- **046** – Accessory rails (*pair/specify head or seat section/factory installed*)
- **450** – Two function foot control (*substitution for hand control/64-*)
- **451** – Simultaneous two function foot and hand control (*69-*)
- **802** – Paper cutter band
- **810** – Restraint straps assembly
- **856** – 220/240 Volt (*factory installed*)
- **5005F** – Fire retardant foam

Table Specifications

- **Height:** Min. 19" (48.3cm) at seat
Max. 35" (88.9cm) at seat
- **Width:** 30" (76.2cm)
- **Length:** 70" (177.8cm) to 80" (203.2cm) with leg section extended
- **Storage drawer:** 16" D x 19" W x 4.5" H (40.6cm D x 48.3cm W x 11.4cm H)
- **Electrical:** 120 VAC, 60 HZ or 220/240 VAC, 60HZ

Colors

